



## SPECIFICATION =

- dispersion of solid particles of size ranging from 1 μm to 50 μm (aerodynamic diameter)
- Venturi tube flow rate is 16.5 litre/min
- output particle concentration from 0.3 to 40 mg/m³
- flow rate of carrier gas ranging from 12 to 21 litre/min
- accessories: device for cleaning, dehumidification and pressure control of the input air





INVESTMENTS IN EDUCATION DEVELOPMENT

## Small-scale powder disperser TSI - SSPD 3433

TSI Model 3433 Small-Scale Powder Disperser is a compact and portable device for dry dispersion of solid particles in gaseous environment. It employs a turntable with a Venturi tube for dispersion of small (milligrams) amounts of dry powder. A Venturi tube creates sufficient shear forces to deagglomerate entering aerosol particles. The disperser generates high quality and time-stable particle dispersion of a defined size range and concentration.





## TYPICAL APPLICATIONS

- generation of aerosols for toxicological studies or for pharmaceutical industry
- generation of tracer particles for gas flow studies
- dispersion of calibrated particles (PSL) for validation measurements of particle sizers
- dispersion of powder particles for their size measurement
- re-dispersion of particles deposited on membrane filters intended for air quality monitoring
- generation of solid particles for research on deposition in environmental engineering

Ing. Jan Jedelský, Ph.D.

tel: +420.541.143.266 | e-mail: jedelsky@fme.vutbr.cz

Department of Thermodynamics and Environmental Engineering

Faculty of Mechanical Engineering, Brno University of Technology

tel: +420 541 143 280 | e-mail: ottp@fme.vutbr.cz | www.eu.fme.vutbr.